



Application overview


1. Start view of the application - from this view you can add ct and rt struct structures files or exit the application.




2. Adding directory with ct files




3. Loading rt struct structures file




4. Informing the user about the ongoing process of reading ct photos, along with rt struct structures.



5. Exemplary application view with ct image uploaded, with rt struct structures, based on the preset settings.



6. Changing the contrast of the image by moving the mouse and pressing the left button. Changing the right-left changes the width of the window, while moving up-down changes the height of the window. Two images with different window parameters are presented below.



7. Option to save the file to any format specified by the user. The process of saving the file is presented below, and the exported photo in png format is inserted at the very end.

